

## Pull out Test Summary

Pull out test of the glass fiber post with three different cements performed at the University of Groningen (NL)

### Materials and Methods

Patent™ two-piece implants were embedded in epoxy to be fixed in a universal testing machine. Glass fiber posts were cemented with three different cements, table 1, according to the manufacturer's instructions. After curing, the samples were individually fixed in the universal testing machine, the glass fiber posts were pulled out and the maximum pull force was registered, **figure 1** and **2**. Five samples were tested for each cement.



Figure 1: Test set-up in universal testing machine.




Figure 2: Close-up of the embedded implant with the cemented post.

Manufacturer	Brand
3M ESPE	RelyX
Ivoclar	Multilink
Kerr	Nexus 3

Table 1: Cements tested.

### Results

The results are presented in diagram 1.

**RelyX** demonstrated a pull-out force of  $1391 \pm 28$  N, **Multilink**  $796 \pm 22$  N and **Nexus 3**  $507 \pm 25$  N.

### Conclusion

All three tested cements demonstrated a pull-out resistance at a level that is adequate for clinical use. **RelyX** demonstrated the highest pull-out resistance.

### Pull out Test




Diagram 1: Pull-out force for the different cements.